


# 10. Infrared Spectroscopy

## A. Background and Introduction

The purpose of this laboratory experiment is for you to gain hands-on experience in the operation of an infrared (IR) spectrometer and interpretation of IR spectra to elucidate structural details. This experiment will supplement the theoretical aspects of IR spectroscopy discussed in lecture. The information provided below was written to be concise and is meant to be a basic summary of IR spectroscopy. You should refer to your lecture notes and/or textbook for additional details. It is imperative that you develop a firm foundation of IR spectroscopy early in the semester, as we will frequently use IR spectroscopy to characterize the products prepared in the organic laboratory.

Organic chemists typically use infrared (IR) spectroscopy to identify the functional groups present in a molecule. IR spectroscopy takes advantage of light energy in the infrared region of the electromagnetic spectrum. When an organic molecule is irradiated with IR light, the light energy is absorbed if the frequency of light is equal to the vibrational energy associated with a particular bond in a molecule. Covalent bonds are not static, but are rather like springs with atoms attached. These springs (bonds) can stretch (2 atoms separated by a bond) and can bend (3 atoms separated by two bonds) as shown in figure 1. The stretching and bending frequencies for a particular molecule are quantized, thus they occur at discrete energy levels. This fact results in every distinct molecule having its own individual IR spectrum that can be distinguished from that of another molecule.



**Figure 1. Stretching and Bending Vibrational Modes**

An IR spectrum is a plot of the percentage of light transmitted (%T) versus the frequency in wavenumbers ( $\tilde{\nu}$ ,  $\text{cm}^{-1}$ ). Figure 2 shows the IR spectrum of ethanol ( $\text{CH}_3\text{CH}_2\text{OH}$ ). Notice the typical IR frequency scale ranging from 400 to 4000  $\text{cm}^{-1}$ . When observing an IR spectrum, it is useful to visualize an imaginary line at 1500  $\text{cm}^{-1}$ . To the right of this line is the “fingerprint region”, which is typically difficult to analyze. To the left of this line is the “functional group region”, which is the most useful for our purposes. The leftmost absorption (circled) represents the O-H bond stretch. This absorption occurs at a relatively high wavenumber (3350  $\text{cm}^{-1}$ ). Analysis of this absorption reveals a 10% transmittance of light, which means 90% of the light energy was absorbed by the molecule at that particular frequency.





Figure 2. Sample IR Spectrum (Ethanol)


Further application of the spring analogy reveals the relationship between the frequency of IR light absorption, bond strength, and atom mass. Hooke's law describes the motion of a vibrating spring where  $k$  is a constant,  $f$  is the force constant, and  $m$  is the mass.

$$\tilde{\nu} = k \sqrt{\frac{f}{m}} \quad (\text{Hooke's Law})$$


Because the force constant,  $f$ , represents the strength of the spring (bond), a stronger bond gives rise to a larger  $\tilde{\nu}$ , and thus represents a higher energy absorption. For this reason, triple bonds absorb at a higher wavenumber than double bonds, which absorb at a higher wavenumber than single bonds.



The mass,  $m$ , of the connected atoms also affects the wavenumber of the absorption. Lighter atoms increase the value of  $\tilde{\nu}$ , which corresponds to a higher energy absorption. For example, the C-H absorption occurs at a higher energy (wavenumber) than does the C-O absorption because hydrogen is a much lighter atom than oxygen.



It is also important to understand that in order for a bond to absorb IR light, there must be a change in the bond's dipole moment during the molecular vibration. For example, 2,3-dimethyl-2-butene has an **IR inactive** C=C bond stretch due to molecular symmetry. During the C=C bond stretch there is no change in the bond's dipole moment.



## Use of the Spectrometer

Our laboratories use an Attenuated Total Reflectance Fourier Transform IR (ATR FT-IR) spectrometer. This type of instrument requires little sample preparation and can easily record the IR spectrum of solids, liquids, and oils. Your TA will show you how to properly use this instrument to acquire an IR spectrum. Although little sample preparation is required, it is important that your sample be free of solvent and other impurities. If your compound of interest contains solvent, and a spectrum is recorded, your spectrum will be contaminated with IR peaks corresponding to the solvent.

## Interpretation of IR Spectra

When given an IR spectrum to interpret, you can use a chart such as the one shown in figure 3 to help you pick out important absorptions corresponding to various functional groups in the molecule. Additionally, the absence of IR absorptions in a particular region gives you important information indicating the absence of certain functional groups.




Figure 3. Some Common IR Regions

### A Detailed Look at Some Important Absorptions








#### 1. Hydroxyl O-H Stretch and Amine N-H Stretch

Structural Unit	Wavenumber, $\text{cm}^{-1}$	Special Features
	3200-3600 (s, br)	
	3300-3500 (m)	$\text{R}_2\text{N-H}$ = one IR stretch $\text{RNH}_2$ = two IR stretches

The hydroxyl O-H stretch is typically broad and strong and appears in the 3200 – 3600  $\text{cm}^{-1}$  range. In the case of a carboxylic acid, the O-H stretch is typically very broad and appears in the 2500 – 3300  $\text{cm}^{-1}$  range, often overlapping the C-H bond stretches. A carboxylic acid will also have a C=O stretch in the 1700 – 1735  $\text{cm}^{-1}$  range.

Amine N-H stretches typically appear in the 3300-3500  $\text{cm}^{-1}$  range and are much more narrow and less intense than an -OH stretch. Primary amines have two amine hydrogen and exhibit two N-H stretches. Secondary amines have only one amine hydrogen and exhibit only one N-H stretch. Tertiary amines have no absorption in this region because there are no hydrogen atoms on the amine nitrogen. In the case of an amide, the NH stretches appear in the same 3300-3500  $\text{cm}^{-1}$  range. Additionally, look for the amide C=O stretch in the 1650 – 1700  $\text{cm}^{-1}$  range.


## 2. Carbonyl C=O Stretch

Structural Unit	Wavenumber, $\text{cm}^{-1}$	Special Features
	<b>1650-1850 (s)</b>	Variable depending on the carbonyl functionality (see below)
	<b>1750-1850</b>	
	<b>1700-1750</b>	Also look for strong $\text{Csp}^2\text{-O}$ stretch between 1200 and 1300 $\text{cm}^{-1}$
	<b>1700-1735</b>	Also look for a strong O-H stretch between 2500 and 3300 $\text{cm}^{-1}$
	<b>1720-1740</b>	Also look for aldehyde C-H stretches at $\sim 2720$ and $\sim 2820$ $\text{cm}^{-1}$
	<b>1680-1750</b>	Generally around 1720 $\text{cm}^{-1}$ ; Decreased when in conjugation
	<b>1650-1700</b>	Also look for N-H stretch between 3300 and 3500 $\text{cm}^{-1}$

The carbonyl region, 1650-1850  $\text{cm}^{-1}$ , is one of the most important regions of the spectrum. The table above lists the most common carbonyl containing functional groups and the relative placement of the C=O stretch in each. Use these regions as a guide, but be aware that these cutoffs are not strict and the C=O stretch for one of these functional group may lie outside of the listed range.

Take note of the special features listed for each of the carbonyl containing functional groups, which can be used to help you confirm its identity. For example, it is very difficult to distinguish the C=O stretch of aldehydes and ketones, however, in the case of an aldehyde, there will also be two aldehyde C-H stretches at  $\sim 2720$  and  $\sim 2820$   $\text{cm}^{-1}$ .

When a carbonyl group is **conjugated** with a double bond or benzene ring, the C=O stretch will be decreased by 20-30  $\text{cm}^{-1}$  (figure 4).



**Figure 4. Effect of Conjugation on Carbonyl Stretching Frequency**

### 3. Alkyne and Nitrile stretches

Structural Unit	Wavenumber, $\text{cm}^{-1}$	Special Features
	<b>2200-2300 (s)</b>	
	<b>~2150 (v)</b>	For terminal alkynes, look for C-H stretch at $\sim 3300 \text{ cm}^{-1}$

The  $\text{C}\equiv\text{C}$  and  $\text{C}\equiv\text{N}$  stretches can be very difficult to distinguish because they both appear around  $2150 \text{ cm}^{-1}$ . Terminal alkynes can be distinguished by the additional C-H stretch at  $\sim 3300 \text{ cm}^{-1}$ .

### 4. Alkene and Aromatic C=C Stretch

Structural Unit	Wavenumber, $\text{cm}^{-1}$	Special Features
	<b>1600-1700 (v)</b>	
	<b>1450-1600 (v)</b>	Look for 2 or 3 peaks in this region generally at $\sim 1600$ , $\sim 1500$ , and $<1500 \text{ cm}^{-1}$

The  $\text{C}=\text{C}$  stretch varies in intensity, but is often rather weak. Typically  $\text{C}=\text{C}$  bond stretches occur between  $1600$  and  $1700 \text{ cm}^{-1}$ , however, if the  $\text{C}=\text{C}$  is in conjugation with a  $\text{C}=\text{O}$ , the  $\text{C}=\text{C}$  stretching frequency will be lowered by  $20\text{-}30 \text{ cm}^{-1}$ .

Aromatic compounds show several  $\text{C}=\text{C}$  absorptions in the  $1450 - 1600 \text{ cm}^{-1}$  region. Typically you will see absorptions at  $\sim 1600$ ,  $\sim 1500$ , and  $\sim 1430 \text{ cm}^{-1}$ . Sometimes you may only be able to see 2 of these absorptions.


### 5. C-H Stretching and Bending

Structural Unit	Wavenumber, $\text{cm}^{-1}$	Special Features
	<b>~3300 (s)</b>	
	<b>2850-2960 (m)</b>	Look for $\text{C}_{\text{sp}^3}\text{-H}$ stretches just below $3000 \text{ cm}^{-1}$ H-C-H bending just above $1400 \text{ cm}^{-1}$ - $\text{CH}_3$ bending just below $1400 \text{ cm}^{-1}$
	<b>3000-3100 (m)</b>	Look for $\text{C}_{\text{sp}^2}\text{-H}$ stretches just above $3000 \text{ cm}^{-1}$

The frequency of the C-H stretch has a large dependency on the hybridization of the carbon bonded to the hydrogen. In the case of terminal alkynes, the carbon is  $sp$ -hybridized and the C-H stretch occurs at  $\sim 3300 \text{ cm}^{-1}$ . The C-H stretch for  $sp^2$ -hybridized carbon atoms occurs at a slightly lower frequency ( $3000 - 3100 \text{ cm}^{-1}$ ). Finally, the C-H stretch for  $sp^3$ -hybridized carbon atoms occurs in the  $2850 - 2960 \text{ cm}^{-1}$  range.  $3000 \text{ cm}^{-1}$  is a convenient dividing line to know. Just above  $3000 \text{ cm}^{-1}$  are the  $\text{C}_{\text{sp}^2}\text{-H}$  stretches while just below  $3000 \text{ cm}^{-1}$  are the  $\text{C}_{\text{sp}^3}\text{-H}$  stretches.

In the case of saturated hydrocarbons, C-H bending vibrations can be observed in the  $1375 - 1475 \text{ cm}^{-1}$  range. Methyl ( $-\text{CH}_3$ ) bending is observed just below  $1400 \text{ cm}^{-1}$  while methylene ( $-\text{CH}_2-$ ) and methine ( $\text{R}_3\text{C-H}$ ) bending is observed just above  $1400 \text{ cm}^{-1}$ .

## 6. C-O Bond Stretch

Structural Unit	Wavenumber, $\text{cm}^{-1}$	Special Features
	<b><math>\text{Csp}^3\text{-O}</math>: 1000-1100 (m)</b> <b><math>\text{Csp}^2\text{-O}</math>: 1200-1300 (s)</b>	Common for alcohols and ethers. Common for esters.

Despite being in the fingerprint region, the C-O single bond stretch is exceptionally useful. The C-O stretch is of medium to strong intensity and appears in the  $1000 - 1300 \text{ cm}^{-1}$  region. The hybridization of the carbon bonded to the oxygen has a large effect on the C-O stretching frequency. When the carbon is  $\text{sp}^3$ -hybridized such as that in alcohols and alkyl ethers, the C-O stretch absorbs in the  $1000 - 1100 \text{ cm}^{-1}$  range. When the carbon is  $\text{sp}^2$ -hybridized, however, the C-O stretch absorbs in the  $1200 - 1300 \text{ cm}^{-1}$  range.

### Analysis of a Sample Spectrum

In the analysis of IR spectra, it is important to get a feel for the general location of the major functional group absorptions. Then, when confronted with an IR spectrum, such as the one shown below, you will be able to examine it quickly and immediately pick out important absorptions corresponding to various functional groups. This information will allow you to both confirm an expected organic structure and to help elucidate the structure of an unknown molecule.








Figure 5. IR Spectrum of Methyl 3-Hydroxybutanoate

## B. Experimental Procedure

Your TA will assign you two different unknowns. Record the IR spectrum of each unknown. Go through each IR spectrum and identify as many functional groups as possible. In your laboratory notebook make a list of characteristic IR absorptions for each compound along with the associated functional group. Draw two potential structures for each unknown that account for all of the indicated functional groups. Based on the limited data you have from the IR spectrum, it is highly unlikely that you will be able to determine the exact structure of your unknown.

### C. Pre-Lab Questions

1. Carbon dioxide (CO<sub>2</sub>) has two types of bond stretches: symmetrical and unsymmetrical. Classify each one of these stretches as IR active or IR inactive and explain your choice.
2. How could you use IR spectroscopy to differentiate between the two isomers: 1-butyne and 2-butyne?
3. Which bond is stronger: the C=O of an ester (1735 cm<sup>-1</sup>) or the C=O bond of a ketone (1715 cm<sup>-1</sup>). Explain your answer.
4. For each compound below, approximate the most important IR absorptions that you would expect to see.



5. For each group of IR frequencies listed below, suggest the functional group that is present.
  - a) 1734, 1250, 1080 cm<sup>-1</sup>
  - b) 3400 (broad), 1050 cm<sup>-1</sup>
  - c) 3050, 1650 cm<sup>-1</sup>

### D. Post-Lab Questions

1. For each of your IR unknowns, make a list of the most important absorbance bands and along with the identity of each absorbance (i.e. 1650 cm<sup>-1</sup>, C=C stretch).
2. Propose two different potential structures for each of your unknowns that are consistent with the IR data that you collected.
3. Explain how IR spectroscopy could be used to monitor the hydroboration-oxidation reaction of 1-hexene to give 1-hexanol.

### E. Acknowledgements

The IR spectra used throughout this experiment were taken from SDBSWeb: <http://sdb.sdb.aist.go.jp> (National Institute of Advanced Industrial Science and Technology, 2016).